





Monteris Medical Solutions Portfolio



The Monteris Medical portfolio of solutions includes products to support patient stabilization and transport from OR to MRI (AtamA™), multiple ways to manage laser probe trajectory alignment (AXiiiS® and Monteris® Mini-Bolt), therapy planning and analysis software (NeuroBlate® M*Vision™) and a versatile set of tools for robotic laser thermotherapy (NeuroBlate® Robotic Probe Driver, SideFire™ and FullFire™ laser probes). For more details on each product, please see individual inserts.

NeuroBlate® System Layout



Not every neurosurgical lesion is small, easy to access and symmetrical

Corporate Headquarters
16305 36th Avenue North Suite 200,
Plymouth, MN 55446
Phone: (763) 253-4710
Fax: (763) 762-2720
marketing@monteris.com

Canadian Office
Unit 1B – 25 Scurfield Blvd.,
Winnipeg, MB R3Y 1G4 Canada
Phone: +1 (204) 272-2220
Fax: +1 (204) 272-2219
marketing@monteris.com



Helping you manage the full range of patient brain abnormalities

Aside from wishing patients had no brain tumors or lesions, it would be great if the ones they did present with were small, symmetrical and easily accessible. But because that's not the case, Monteris is there to provide you a robust set of tools to address various sizes, densities, shapes and locations of your patient neurosurgical lesions.

With high capacity and directional lasing per choice, NeuroBlate can effectively cover wide range lesions. For a symmetrical lesion, using NeuroBlate® FullFire Diffused Tip Probe can be efficient. Monteris measures performance as ablated lesion volume per minute, which averages 1.64 cc per min.* For a 2 cm diameter or less lesion, it takes as little as 2 minutes. However, NeuroBlate technology provides you versatility, in addition to efficiency, enabling you to address difficult and complex patient circumstances. It's about time...that you have the tools you need for the quickly ablated abnormalities, coupled with the flexible toolset of the full Monteris Medical solutions portfolio.



*Monteris data on file.

Product Support Services

In addition to product innovation and reliability, Monteris Medical backs your NeuroBlate® program with these additional services:

- Reimbursement Support – Dictation tools, coding guides, and claims assistance through a third party resource are available for your NeuroBlate® program,
- NeuroBlate® Outreach Services – The NeuroBlate® Outreach program offers many resources and templates to help you reach out effectively to patients, referring physicians and the broader community,
- NeuroBlate® University – A series of educational events, resources and web content for care team use before, during and after the NeuroBlate® procedure, and
- NeuroBlate® Network – Enrollment in the Network, which is a community of NeuroBlate® surgeons and their facilities to interact with one another, sharing interesting cases, best practices, and providing ongoing mentoring support.

Simple to Complex...Innovation for Multiple Tumor Types


LESION TYPES							
Not all tumors are created equal							
	2 cm or Less	Greater than 2 cm, Symmetrical	Greater than 2 cm, Asymmetrical	Difficult to Access	Temporal Lobe	Multiple	Posterior Fossa

PROBE CHOICES							
NeuroBlate®							
FullFire™	•	•		•	•	•	•
SideFire™		•	•	•	•	•	•

Shown above are suggested solutions from Monteris Medical, depending on the characteristics of the brain lesion you are ablating.

NeuroBlate® FullFire™ and FullFire™ Select Diffusing Tip Laser Probe

NeuroBlate® SideFire™ and SideFire™ Select Directional Laser Probe



NOTE: It is for the physician to determine which product to use for which patient condition. This tool is provided as a general guide for convenience of product appropriateness.

