

# NeuroBlate®

S Y S T E M

## ATAMA SYSTEM

PATIENT TRANSFER &  
HEAD STABILIZATION  
SYSTEM FOR THE MRI



Enables attachment of NeuroBlate  
System hardware\*

Integrates versatile head fixation into a  
patient transfer board for efficient patient  
transport from OR to MRI

1.5 or 3.0 T MR conditional



Secured patient head position ensures temporal and spatial  
imaging accuracy for MRI-guided neuro-interventions

Starburst Adapter enables attachment of image  
guided surgery hardware for navigation


\* NeuroBlate System offers Robotic Laser Thermotherapy

# SPLIT RING HEAD FIXATION


- Removable head fixation ring allows stereotactic frame-like head immobilization.
- Head fixation ring can be rotated and fixed for flexible positioning.
- Multiple positions around the ring and on fixation posts are available to provide flexible pinning options.



Split ring design provides simplified application and quick removal



Posterior section of ring can be used alone




Extended front ring section accommodates larger heads

## DISCLOSURES:

The NeuroBlate® System is intended for ablating intracranial soft tissue, including brain structures. Patients must be able to undergo MRI exposure and be surgical candidates. The technology is not appropriate for every lesion type and location. It may be difficult to use the technology on certain large or irregularly shaped lesions.

Possible adverse events include, but are not limited to, compromised device function, hematoma, embolic events, edema, bleeding, unintended major tissue damage and permanent neurological deficits. Prior to using these devices, please review the Instructions for Use for a complete listing of indications, contraindications, warnings, precautions and potential adverse events. For full prescribing information, please visit [monteris.com](http://monteris.com).



## ATAMA™ SYSTEM CONFIGURATIONS:

ATB002: Siemens (Espree or Avanto), IMRIS, and GE MRI

ATB002-AS: Siemens Aera or Skyra

ATB002-AC: Philips Achieva

ATB002-IG: Philips Ingenia

## ATAMA SYSTEM INCLUDES:

- Patient Transfer Board
- Rolling Storage Stand
- ATF102 & ATF103: Head Fixation Ring Assembly (two & three piece) with Tray (also sold separately)
- ASA001: Starburst Adapter for Trajectory Devices with a Standard Starburst Interface

Contact Monteris Medical Corporation for more information. Not available for sale outside the U.S. or Canada.

### U.S. Office

14755 27th Ave. N., Suite C  
Plymouth, MN 55447, USA  
+1.866.799.7655

[monteris.com](http://monteris.com)

Rx Only

### Canadian Office

Unit 1B-25 Scurfield Blvd.  
Winnipeg, MB R3Y1G4, Canada  
+1.204.272.2220

