


MONTERIS MINI-BOLT

CRANIAL BOLT SYSTEM WITH
ROBOTIC INTERFACE



Rigid skull fixation for neurosurgical devices

Solid titanium construction

Available in 2.2 mm and 3.3 mm
inner diameter

4.5 mm diameter hole

Allows a direct interface to the NeuroBlate® Robotic
Probe Driver for precise laser probe control and
laser delivery

Accessories and adapters allow on-trajectory
placement using a wide array of stereotactic frames
or frameless navigation and robotic systems

1.5 or 3.0 T MR Conditional




SOLUTIONS FOR THE FULL SPECTRUM OF STEREOTACTIC DELIVERY OPTIONS

The only cranial bolt system that enables a robotic interface for protected and precise laser delivery.

Efficiency – Trimmed-body-and-tapered-thread design provides time-saving deployment and compact spacing for multiple trajectories

Confidence – Solid titanium construction sustains stereotactic accuracy for MRI-guided neurosurgical interventions

Adaptation – Array of accessories supports on-trajectory placement using the most common stereotactic frames, articulated arms and robot



MONTERIS MINI-BOLT ACCESSORIES

Contact your Monteris Medical representative to order individual kit components.

Monteris® Mini-Bolt	
CMB033	Monteris Mini-Bolt for 3.3 mm Device
CMB022	Monteris Mini-Bolt for 2.2 mm Device
Monteris® Mini-Bolt Accessory Kits	
Alignment Mandrel, Insert Adapter for Alignment Mandrel, 4.5 mm Clearance Mandrel, 4.5 mm Depth Stop, 4.5 mm Insert Adapter, T-Handle with Mini-Bolt Stereotactic Driver, Mini-Bolt Removal Tool for MRI, Sterilization Tray with Dividers	
CMB033-AA	3.3 Monteris Mini-Bolt Accessory Kit for 3.3 mm Mini-Bolt
CMB022-AA	2.2 Monteris Mini-Bolt Accessory Kit for 2.2 mm Mini-Bolt
CMB033-UP	Upgrade to 3.3 mm Mini-Bolt Accessories (3.3 mm Alignment Mandrel, Insert Adapter and Stereotactic Driver only)
CMB022-UP	Upgrade to 2.2 mm Mini-Bolt Accessories (2.2 mm Alignment Mandrel, Insert Adapter and Stereotactic Driver only)
Additional Stereotactic Device Adapters and Accessories	
CMB-CW	Accessory Host Adapter for CRW Stereotactic Frame
CMB-LK	Accessory Host Adapter for Leksell Stereotactic Frame
CMB-RR	Accessory Host Adapter for Rosa Robot
CMB-MA	Accessory Host Adapter for Autoguide System
CMB022-SD	6 mm Specialized Driver for 2.2. Mini-Bolt with Insert Adapter
CMB-SD	6 mm Specialized Driver for 3.3 Mini-Bolt with Insert Adapter
CMB-34	Insert Adapter for 3.4 mm ID
CMB-TA	Trumpet Insert Adapter 4.5 mm ID
CMB-TR	Mini-Bolt Sterilization Tray with Dividers

DISCLOSURES:

Prior to using these devices, please review the Instructions for Use for a complete listing of indications, contraindications, warnings, precautions, and potential adverse events. For full prescribing information, please visit monteris.com.

Contact Monteris Medical Corporation for more information. Not available for sale outside the U.S. or Canada. Monteris Mini-Bolt and AXiiiS-CMB are identical products.

U.S. Office

14755 27th Ave. N., Suite C
Plymouth, MN 55447, USA
+1.866.799.7655

Canadian Office

Unit 1B-25 Scurfield Blvd.
Winnipeg, MB R3Y1G4, Canada
+1.204.272.2220

monteris.com

Rx Only

©2020 Monteris. All Rights Reserved. TM and [®] denote Trademarks and Registered Trademarks of Monteris Medical Corporation



AA10065 Rev G 7/20