



MINI-BOLT

CRANIAL BOLT WITH ROBOTIC INTERFACE

MINIMALLY INVASIVE

- Available in 2.2 mm and 3.3 mm inner diameter
- 4.5 mm diameter hole

PRECISE

- Rigid skull fixation for neurosurgical devices
- Precise laser probe placement on-trajectory
- Direct interface with NeuroBlate Robotic Probe Driver for added control

EFFICIENT WORKFLOW

- Compatible with most stereotactic frames and navigation systems
- Supplied sterile in a dual-barrier pouch
- 1.5T and 3.0T MR Conditional

SOLUTIONS FOR A RANGE OF STEREOTACTIC DELIVERY OPTIONS

Mini-Bolt and Accessories

Contact your Monteris Medical representative to order individual kit components.

Mini-Bolt

CMB033-S	Mini-Bolt for 3.3 mm Device - Sterile
CMB022-S	Mini-Bolt for 2.2 mm Device - Sterile
CMB033	Mini-Bolt for 3.3 mm Device
CMB022	Mini-Bolt for 2.2 mm Device

Mini-Bolt Accessory Kits

Alignment Mandrel, Insert Adapter for Alignment Mandrel, 4.5 mm Clearance Mandrel, 4.5 mm Depth Stop, 4.5 mm Insert Adapter, T-Handle with Mini-Bolt Stereotactic Driver, Mini-Bolt Removal Tool for MRI, Sterilization Tray with Dividers

CMB033-AA	3.3 Mini-Bolt Accessory Kit for 3.3 mm Mini-Bolt
CMB022-AA	2.2 Mini-Bolt Accessory Kit for 2.2 mm Mini-Bolt
CMB033-UP	Upgrade to 3.3 mm Mini-Bolt Accessories (3.3 mm Alignment Mandrel, Insert Adapter and Stereotactic Driver only)
CMB022-UP	Upgrade to 2.2 mm Mini-Bolt Accessories (2.2 mm Alignment Mandrel, Insert Adapter and Stereotactic Driver only)

Additional Stereotactic Device Adapters and Accessories

CMB-CW	Accessory Host Adapter for CRW Stereotactic Frame
CMB-LK	Accessory Host Adapter for Leksell Stereotactic Frame
CMB-RR	Accessory Host Adapter for Rosa Robot
CMB-MA	Accessory Host Adapter for Autoguide System
CMB022-SD	6 mm Specialized Driver for 2.2 Mini-Bolt with Insert Adapter
CMB-SD	6 mm Specialized Driver for 3.3 Mini-Bolt with Insert Adapter
CMB-34	Insert Adapter for 3.4 mm ID
CMB-TA	Trumpet Insert Adapter 4.5 mm ID
CMB-TR	Mini-Bolt Sterilization Tray with Dividers

Disclosures:

MiniBolt Cranial Fixation Bolt, also known as AXiiiS CMB. The Mini-Bolt is a disposable device intended to provide placement and skull fixation of neurosurgical instruments or devices with an outer diameter (OD) up to 3.3 mm or 2.2 mm.

All brain surgeries carry risk. Possible adverse events include, but are not limited to, tissue injury, infection, and cerebral spinal fluid (CSF) leak. Prior to using these devices, please review the Instructions for Use for a complete listing of indications, contraindications, warnings, precautions and potential adverse events. For full prescribing information, please visit monteris.com.

Rx Only.

U.S. Office

131 Cheshire Lane, Suite 100
Minnetonka, MN 55305
+1.866.799.7655

Canadian Office

Unit 1B-25 Scurfield Blvd.
Winnipeg, MB R3Y1G4, Canada
+1.204.272.2220



monteris.com

Rx Only