



NeuroBlate® System and Atama Instructions for Use Compatible MRI Systems Supplement

80575 Rev B – Sept. 2024

Contents

1. Overview	1
1.1. Warning and Identification Labels	2
2. NeuroBlate Compatible MRI Systems.....	2

Patent Notice under 35 U.S.C. §287(a): Monteris Medical products and associated processes are subject to the United States patents and pending applications.

NeuroBlate, Monteris Medical, SideFire, FullFire, Atama, AutoLITT, Monteris Mini-Bolt, the M Logo, the N Logo, the A Logo and Fusion-S Pro are trademarks of Monteris Medical. ©2023 Monteris. All Rights Reserved. ™ and ® denote Trademarks and Registered Trademarks of Monteris Medical Corporation.

1. Overview





This document describes the MRI Systems that are compatible with the NeuroBlate System.

This document supplements information for the following NeuroBlate System IFUs:

- > 80577 – NeuroBlate System v3.18 SW
- > 80554 – NeuroBlate System v3.17 SW
- > 80512 – NeuroBlate System v3.16 SW
- > 80327 – Atama System

1.1. WARNING AND IDENTIFICATION LABELS

Symbols displayed on Monteris products or in their documentation are:

	<p>MR Unsafe - item is NOT MRI compatible and is known to pose a hazard in MR environments. This equipment should not be taken into the MRI room within the 5 Gauss perimeter line.</p>
 <p>1.5T / 3.0T</p>	<p>MR Conditional - the item has been demonstrated to pose no known hazards in a 1.5T or 3.0T MR environments as described in Section 4.4 of the Instructions for Use.</p>
	<p>MR Safe - the item poses NO known hazards in ALL MR environments.</p>
	<p>Refer to instructions</p>

2. NeuroBlate Compatible MRI Systems

The NeuroBlate System and Atama is compatible with only the following 1.5T or 3.0T MRI systems.



1.5T / 3.0T

- > Siemens - minimum software versions VB17, VB19, VD13, VE11, (VXA20, VXA30, VXA50, VXA51, VXA60 and VXA61 with Access-I enabled systems)
 - Symphony (1.5T)
 - Avanto (1.5T)
 - Trio (3.0T)
 - Espree (1.5T)
 - Verio (3.0T)
 - Aera (1.5T)
 - Skyra (3.0T)
 - Prisma (3.0T)
 - Sola (1.5T)
 - Vida (3.0T)
- > IMRIS NEURO II (1.5T) or NEURO III (3.0T) - Siemens MRI software versions listed above
- > GE - minimum software version 15.X to 30.X
 - Signa HDx/HDxT (1.5T & 3.0T)



U.S. Office

131 Cheshire Lane, Suite 100
Minnetonka, MN 55305, USA
+1.866.799.7655

Canadian Office

Unit 1B-25 Scurfield Blvd.
Winnipeg, MB R3Y1G4, Canada
+1.204.272.2220



- Optima MR450/450W (1.5T)
- Discovery MR750/750W (3.0T)
- Signa Artist (1.5T)
- Signa Architect (3.0T)
- Signa Premier (3.0T)

- > Philips – software version R3, R4, R5
 - Achieva (1.5T & 3.0T)
 - Ingenia (1.5T and 3.0T)

- > Intera (1.5T) software version v3



U.S. Office

131 Cheshire Lane, Suite 100
Minnetonka, MN 55305, USA
+1.866.799.7655

Canadian Office

Unit 1B-25 Scurfield Blvd.
Winnipeg, MB R3Y1G4, Canada
+1.204.272.2220

AA80575 Rev B 9/24

